





The EZ Rom Elite adjustable knee brace allows the user to adjust the length of the brace for custom fitting. This product is ideal for knee ligament damage, knee surgeries, and instability of the knee. The EZ Rom Elite is designed for post-operative fixation, allowing the doctor to lock the patient into extension and control range of motion as prescribed.

Post-Op Adjustable Knee Brace

Features and Benefits

Telescoping and malleable hinge bars extendable from 16.5" – 27" (42cm – 68.5cm) Hinge range of motion adjustment: Flexion (0° to 120°) and extension (0° to 90°). Quick lock key offers ability to control and lock out the hinge motion to 0°, 15° and 30° directly. Light weight and breathable material. Easy application with hook and buckle closure. Latex free. Universal size, design fits either left or right leg.

Indications for Use

ACL, PCL, MCL, LCL repairs/injuries Tibial plateau fractures Osteochondral repairs Meniscal repairs

Patella tendon repairs Sprains/strains of the knee HTOs (high Tibial osteotomy)



Safety Lock



Hinge with quick lock feature

Order information Part Number: EZRE001

Size: Universal

How to order Email: orders@manamed.net, Fax: 888-508-0712, Call: 888-508-0712







Extension

To control Extension, pull the extension tab out and move to desired setting between 0° and 90° and release tab to lock into position.

Flexion

To control Flexion, pull the Flexion tab out and move to desired setting between 0° and 120° and release tab to lock into position.

Quick Lock

The hinge can be locked out at $0^{\circ}, 15^{\circ}$, 30° by sliding the yellow tab into the locked position.